

## **FREE T3**

**REF A13422**

---

**Intended Use** The Access Free T3 assay is a paramagnetic particle, chemiluminescent immunoassay for the quantitative determination of free triiodothyronine levels in human serum and plasma using the Access Immunoassay Systems.

---

**Summary and Explanation** The hypothalamic pituitary-thyroid axis controls thyroid hormone synthesis, release, and action. Thyrotropin-releasing hormone (TRH) secreted from the hypothalamus stimulates the synthesis and release of thyrotropin or thyroid-stimulating hormone (hTSH). hTSH, in turn, stimulates the synthesis, storage, secretion, and metabolism of thyroxine (T4) and triiodothyronine (T3).

T3 is the major biologically active thyroid hormone. Of the circulating T3, about 80% is formed from peripheral deiodination of thyroxine and 20% is secreted directly from the thyroid gland. The T4 and T3 hormones are transported in the circulation bound to thyroxine binding globulin (TBG), thyroxine binding pre-Albumin (TPBA) and albumin. About 0.2 to 0.4% of the circulatory total T3 is in equilibria as unbound or free, in contrast to about 0.03% of the total T4. In most individuals, the free fractions of these hormones correlate with the functional thyroid state.<sup>1,2</sup>

Free T4 and T3 regulate normal growth and development by maintaining body temperature and stimulating calorogenesis. In addition, free T4 and free T3 affect all aspects of carbohydrate metabolism as well as certain areas of lipid and vitamin metabolism. Fetal and neonatal development also require thyroid hormones.<sup>1,2</sup>

With normal levels of thyroid binding proteins, free T3 levels correlate with total T3. Measuring free T3 is useful when altered levels of total T3 occur due to changes in thyroid hormone binding proteins, especially in cases with altered TBG or low albumin concentrations. Free T3 is elevated alone (T3 toxicosis) in about 5% of hyperthyroids.

Non-thyrometabolic disorders may cause abnormal free T3 levels. Determination of thyroid status in patients with non-thyroidal illness (NTI) should be interpreted with caution.<sup>2,3</sup> For example, anticonvulsant drug therapy (particularly phenytoin) may result in decreased free T3 levels due to an increased hepatic metabolism, and secondarily to displacement of hormone from binding sites.<sup>2,4,5</sup> Anti-inflammatory drugs such as salicylate and phenylbutazone also compete for hormone binding sites, but their effect on free T3 levels has not been clearly defined.<sup>2,6</sup> Patients on heparin therapy may have elevated free T3 levels due to release of non esterified fatty acids (NEFA), which can alter the relationship between free and bound hormones.<sup>5</sup>

---

**Principles of the Procedure** The Access Free T3 assay is a competitive binding immunoenzymatic assay. A sample is added to a reaction vessel with an anti-T3 monoclonal antibody conjugated to alkaline phosphatase. During the incubation, free T3 in the sample reacts with the anti-T3 antibody. Particles coated with streptavidin and biotinylated T3 analog are then added to the mixture. Unoccupied binding sites on the anti-T3 antibody are bridged to the particle through the T3 analog. After incubation in a reaction vessel, materials bound to the solid phase are held in a magnetic field while unbound materials are washed away. Then, the chemiluminescent substrate Lumi-Phos\* 530 is added to the vessel and light generated by the reaction is measured with a luminometer. The light production is inversely proportional to the concentration of free T3 in

the sample. The amount of analyte in the sample is determined from a stored, multi-point calibration curve.

## Product Information

### Access Free T3 Reagent Pack

Cat. No. A13422: 100 determinations, 2 packs, 50 tests/pack

- Provided ready to use.
- Store upright and refrigerate at 2 to 10°C.
- Refrigerate at 2 to 10°C for a minimum of two hours before use on the instrument.
- Stable until the expiration date stated on the label when stored at 2 to 10°C.
- Stable at 2 to 10°C for 28 days after initial use.
- Signs of possible deterioration are a broken elastomeric layer on the pack or control values out of range.
- If the reagent pack is damaged (i.e., broken elastomer), discard the pack.
- All antisera are polyclonal unless otherwise indicated.

<b>R1a:</b>	Dynabeads** paramagnetic particles coated with streptavidin in a TRIS buffer with protein (aves), surfactant, < 0.1% NaN <sub>3</sub> and 0.1% ProClin*** 300.
<b>R1b:</b>	MES buffer and 0.1% ProClin 300.
<b>R1c:</b>	Biotinylated T3 analog in a TRIS buffer with protein (aves), surfactant, < 0.1% NaN <sub>3</sub> and 0.1% ProClin 300.
<b>R1d:</b>	TRIS buffer containing animal protein (goat, bovine, aves), surfactant, < 0.1% NaN <sub>3</sub> , and 0.5% ProClin 300.
<b>R1e:</b>	Monoclonal antibody-alkaline phosphatase conjugate in an ACES buffer with protein (aves), surfactant, < 0.1% NaN <sub>3</sub> and 0.1% ProClin 300.

## Warnings and Precautions

- For *in vitro* diagnostic use.
- Patient samples and blood-derived products may be routinely processed with minimum risk using the procedure described. However, handle these products as potentially infectious according to universal precautions and good clinical laboratory practices, regardless of their origin, treatment, or prior certification. Use an appropriate disinfectant for decontamination. Store and dispose of these materials and their containers in accordance with local regulations and guidelines.
- Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal of liquids, flush with a large volume of water to prevent azide build-up.<sup>7</sup>
- Xi. Irritant: 0.5% ProClin 300.



R 43: May cause sensitization by skin contact.

S 28-37: After contact with skin, wash immediately with plenty of soap and water. Wear suitable gloves.

- The Material Safety Data Sheet (MSDS) is available upon request.

## Specimen Collection and Preparation

1. Serum and plasma (heparin) are the recommended samples.
2. Observe the following recommendations for handling, processing, and storing blood samples:<sup>8</sup>
  - Collect all blood samples observing routine precautions for venipuncture.
  - Allow serum samples to clot completely before centrifugation.
  - Keep tubes stoppered at all times.
  - Within two hours after centrifugation, transfer at least 500 µL of cell-free sample to a storage tube. Tightly stopper the tube immediately.

- Store samples tightly stoppered at room temperature (15 to 30°C) for no longer than eight hours.
  - If the assay will not be completed within eight hours, refrigerate the samples at 2 to 8°C.
  - If the assay will not be completed within 48 hours, or for shipment of samples, freeze at -20°C or colder.
3. Use the following guidelines when preparing specimens:
    - Ensure residual fibrin and cellular matter has been removed prior to analysis.
    - Follow blood collection tube manufacturer's recommendations for centrifugation.
  4. Each laboratory should determine the acceptability of its own blood collection tubes and serum separation products. Variations in these products may exist between manufacturers and, at times, from lot to lot.
  5. Thaw samples no more than three times.

---

**Materials Provided** R1 Access Free T3 Reagent Packs

---

**Materials Required But Not Provided**

1. Access Free T3 Calibrators  
Provided at zero and approximately 1, 2, 5, 10 and 30 pg/mL (1.5, 3.1, 7.7, 15 and 46 pmol/L).  
Cat. No. A13430
2. Quality Control (QC) materials: commercial control material
3. Access Substrate  
Cat. No. 81906
4. **Access, Access 2, SYNCHRON LXi:**  
Access Wash Buffer II, Cat. No. A16792  
**UniCel DxI:**  
UniCel DxI Wash Buffer II, Cat. No. A16793

---

**Procedural Comments**

1. Refer to the appropriate system manuals and/or Help system for a specific description of installation, start-up, principles of operation, system performance characteristics, operating instructions, calibration procedures, operational limitations and precautions, hazards, maintenance, and troubleshooting.
2. Mix contents of new (unpunctured) reagent packs by gently inverting pack several times before loading on the instrument. Do not invert open (punctured) packs.
3. Use fifty-five (55) µL of sample for each determination in addition to the sample container and system dead volumes. Refer to the appropriate system manuals and/or Help system for the minimum sample volume required.
4. The system default unit of measure for sample results is pg/mL. To change sample reporting units to the International System of Units (SI units), pmol/L, refer to the appropriate system manuals and/or Help system. To manually convert concentrations to the International System, multiply pg/mL by multiplication factor 1.536 (T3 m.w. equals 651 daltons).

---

**Procedure** Refer to the appropriate system manuals and/or Help system for information on managing samples, configuring tests, requesting tests, and reviewing test results.

---

**Calibration Details** An active calibration curve is required for all tests. For the Access Free T3 assay, calibration is required every 28 days. Refer to the appropriate system manuals and/or Help system for information on calibration theory, configuring calibrators, calibrator test request entry, and reviewing calibration data.

---

**Quality Control** Quality control materials simulate the characteristics of patient samples and are essential for monitoring the system performance of immunochemical assays. Because samples can be

processed at any time in a “random access” format rather than a “batch” format, quality control materials should be included in each 24-hour time period.<sup>9</sup> Include commercially available quality control materials that cover at least two levels of analyte. More frequent use of controls or the use of additional controls is left to the discretion of the user based on good laboratory practices or laboratory accreditation requirements and applicable laws. Follow manufacturer’s instructions for reconstitution and storage. Each laboratory should establish mean values and acceptable ranges to assure proper performance. Quality control results that do not fall within acceptable ranges may indicate invalid test results. Examine all test results generated since obtaining the last acceptable quality control test point for this analyte. Refer to the appropriate system manuals and/or Help system for information about reviewing quality control results.

---

**Results** Patient test results are determined automatically by the system software using a weighted four parameter logistic curve (4PLC) math model. The amount of analyte in the sample is determined from the measured light production by means of the stored calibration data. Patient test results can be reviewed using the appropriate screen. Refer to the appropriate system manuals and/or Help system for complete instructions on reviewing sample results.

---

**Limitations of the Procedure**

1. Samples can be accurately measured within the analytic range of the lower limit of detection and the highest calibrator value (approximately 0.88–30 pg/mL [1.4–46 pmol/L]).
    - If a sample contains less than the lower limit of detection for the assay, report the results as less than that value (i.e., < 0.88 pg/mL [< 1.4 pmol/L]).
    - If a sample contains more than the stated value of the highest Access Free T3 Calibrator (S5), report the result as greater than that value (i.e., > 30 pg/mL [> 46 pmol/L]). Samples can not be diluted for free T3 determinations.
  2. For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays. Additionally, other heterophile antibodies such as human anti-goat or human anti-triiodothyronine antibodies may be present in patient samples.<sup>10,11</sup> Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.
  3. The Access Free T3 results should be interpreted in light of the total clinical presentation of the patient, including: symptoms, clinical history, data from additional tests, and other appropriate information.
  4. Patient with biotin (vitamin H) treatment may have elevated free T3 levels due to displacement of the biotinylated T3 analog from paramagnetic particles coated with streptavidin.
  5. Serum and plasma values should not be used interchangeably. When monitoring a patient, it is suggested to rebaseline when changing sample types.
- 

**Expected Values**

1. Each laboratory should establish its own reference ranges to assure proper representation of specific populations.
  2. Concentrations of free T3 were measured in sera from 200 apparently healthy, ambulatory subjects. The median value was 3.2 pg/mL (4.9 pmol/L) with a 95% non-parametric range of 2.5–3.9 pg/mL (3.8–6.0 pmol/L).
-

**Specific  
Performance  
Characteristics**

**Methods Comparison**

A comparison of 397 values using the Access Free T3 assay on the Access Immunoassay system and a commercially available immunoassay kit gave the following statistical data using Deming calculations:<sup>12</sup>

n	Range of Observations (pg/mL)	Intercept (pg/mL)	Slope	Correlation Coefficient (r)
397	0.8–17.5	0.0828	0.9464	0.9874

**Imprecision**

This assay exhibits total imprecision  $\leq 12\%$  at concentrations greater than 2.0 pg/mL.

One study, using commercially available human serum based control material generating a total of 20 assays, 2 replicates per day, over 20 days provided the following data, analyzed via analysis of variance (ANOVA).<sup>13,14</sup>

Human Serum Control	Grand Mean (n=40) (pg/mL)	Within Run (%CV)	Between Run (%CV)	Total Imprecision (%CV)
Level 1	1.4	6.6	8.0	10.4
Level 2	2.6	2.6	5.1	5.7
Level 3	9.6	5.1	1.3	5.3

**Analytical Specificity/Interferences**

Samples containing up to 20 mg/dL (342  $\mu\text{mol/L}$ ) bilirubin, 10 mg/dL (119  $\mu\text{mol/L}$ ) conjugated bilirubin, lipemic samples containing the equivalent of 3000 mg/dL (33.87 mmol/L) triolein, or up to 500 mg/dL cholesterol and samples containing up to 500 mg/dL (5.0 g/L) of hemoglobin do not effect the concentration of free T3 assayed.

The following substances were added to serum samples. When tested in the Access Free T3 assay, the observed mean percent changes were as follows:

Substance	Analyte Added (mg/dL)	% Change
Aspirin	75	1.3
Sodium Salicylate	75	13.9
Ibuprophen	750	17.2
Acetaminophen	200	-5.0
Phenylbutazone	7.5	3.9
Thiouracil	5	2.1
Phenytoin	10	-0.1
Furosemide	2	15.9
Carbamazepine	12	12.7
Methimazol	0.4	-0.1
Oleic Acid	283	11.4
Linoleic Acid	280	0
D-biotin	0.001	-7.3

The following table describes the cross-reactivity of the assay with substances that are similar in structure to T3. Following CLSI C45A <sup>15</sup> recommendations the cross-reactants were tested in the absence of serum proteins (i.e. in HEPES pH 7.7 buffer containing 0.1% gelatin). The cross reactivities were calculated as the ratio between the concentration of T3 required to displace 50% of the conjugate (IC50) and the concentration of the cross-reactant required to give the same displacement, calculated on a molarity % basis.

Substance	Cross-Reactivity (%)
Reverse T3	< 0.5
Tetraiodothyroacetic acid	< 0.5
D-thyroxine	< 0.5
L-thyroxine	< 0.1
3, 5 diiodothyronine	3.3
Diiodo-L-tyrosine	< 0.01
Monoiodotyrosine	< 0.01
3-3', 5-Triiodothyroacetic acid (TRIAC)	13.0

#### Analytical Sensitivity

The lowest detectable level of free T3 distinguishable from zero (Access Free T3 Calibrator S0) with 95% confidence is 0.88 pg/mL (1.4 pmol/L). This value is determined by processing a complete six point calibration curve, controls, and 10 replicates of the zero calibrator in multiple assays. The analytical sensitivity value is interpolated from the curve at the point that is two standard deviations from the mean measured zero calibrator signal.

## FREE T3 CALIBRATORS

**REF** A13430

**Intended Use** The Access Free T3 Calibrators are intended to calibrate the Access Free T3 assay for the quantitative determination of free triiodothyronine levels in human serum and plasma using the Access Immunoassay Systems.

**Summary and Explanation** Quantitative assay calibration is the process by which samples with known analyte concentrations (i.e., assay calibrators) are tested like patient samples to measure the response. The mathematical relationship between the measured responses and the known analyte concentrations establishes the calibration curve. This mathematical relationship, or calibration curve, is used to convert RLU (Relative Light Unit) measurements of patient samples to specific quantitative analyte concentrations.

**Traceability** The measurand (analyte) in the Access Free T3 Calibrators is traceable to the manufacturer's working calibrators. Traceability process is based on EN ISO 17511.

The assigned values were established using representative samples from this lot of calibrator and are specific to the assay methodologies of the Access reagents. Values assigned by other methodologies may be different. Such differences, if present, may be caused by inter-method bias.

**Product Information** **Access Free T3 Calibrators**  
**Cat. No. A13430: S0–S5, 2.5 mL/vial**

- Provided ready to use.
- Store upright and refrigerate at 2 to 10°C.
- Mix contents by gently inverting before use. Avoid bubble formation.
- Stable until the expiration date stated on the label when stored at 2 to 10°C.
- Vial is stable at 2 to 10°C for four months after initial use.
- Signs of possible deterioration are control values out of range.
- Refer to calibration card for exact concentrations.

<b>S0:</b>	HEPES buffer, protein (bovine), surfactant, < 0.1% NaN <sub>3</sub> , and 0.5% ProClin*** 300.
<b>S1, S2, S3, S4, S5:</b>	HEPES buffer, protein (bovine), surfactant, T3, < 0.1% NaN <sub>3</sub> , and 0.5% ProClin 300.
<b>Calibration Card:</b>	1

**Warnings and Precautions**

- For *in vitro* diagnostic use.
- Patient samples and blood-derived products may be routinely processed with minimum risk using the procedure described. However, handle these products as potentially infectious according to universal precautions and good clinical laboratory practices, regardless of their origin, treatment, or prior certification. Use an appropriate disinfectant for decontamination. Store and dispose of these materials and their containers in accordance with local regulations and guidelines.

- Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal of liquids, flush with a large volume of water to prevent azide build-up.<sup>7</sup>
- Xi. Irritant: 0.5% ProClin 300.



R 43: May cause sensitization by skin contact.

S 28-37: After contact with skin, wash immediately with plenty of soap and water. Wear suitable gloves.

- The Material Safety Data Sheet (MSDS) is available upon request.

---

**Procedure** Refer to the appropriate system manuals and/or Help system for information on calibration theory, configuring calibrators, calibrator test request entry, and reviewing calibration data.

---

**Calibration Details** The Access Free T3 Calibrators are provided at six levels – zero and approximately 1, 2, 5, 10 and 30 pg/mL. Assay calibration data are valid up to 28 days.  
Calibrators run in duplicate.

---

**Limitations of the Procedure** If there is evidence of microbial contamination or excessive turbidity in a reagent, discard the vial.

---



## References

- 1 Gornall, AG, Luxton, AW, Bhavnani, BR. Endocrine disorders. In Applied Biochemistry of Clinical Disorders. 1986, 305-318. Philadelphia, PA: J. B. Lippincott Co.
- 2 White, GH. Recent advances in routine thyroid function testing. CRC - Critical Reviews in Clinical Laboratory Sciences, 1987, 24: 315-362.
- 3 Spencer, CA. Thyroid status: trends in testing – selective test use cuts cost. Clinical Chemistry News, November, 1989, 15(11): 9-14.
- 4 Liewendahl, K, Majuri, H, Helenius, T. Thyroid function tests in patients on long-term treatment with various anticonvulsant drugs. Clinical Endocrinology, 1978; 8: 187-191.
- 5 Wenzel, KW. Pharmacological interference with in vitro tests of thyroid function. Metabolis. 1981; 30(7): 717-732.
- 6 Wilke, TJ. Estimation of free thyroid hormone concentrations in the clinical laboratory. Clinical Chemistry. 1986; 32(4): 585-592.
- 7 DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard. Available <http://www.cdc.gov/niosh>.
- 8 Approved Guideline – Procedures for the Handling and Processing of Blood Specimens, H18-A3. 2004. Clinical and Laboratory Standards Institute.
- 9 Cembrowski GS, Carey RN. Laboratory quality management: QC  $\Rightarrow$  QA. ASCP Press, Chicago, IL, 1989.
- 10 Kricka L. Interferences in immunoassays – still a threat. Clin Chem 2000; 46: 1037–1038.
- 11 Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613–621.
- 12 Cornbleet, JP, Gochman N. Incorrect least-squares regression coefficients in method-comparison analysis. Clinical Chemistry, 1979, 25(3): 432-438.
- 13 Approved Guidelines - Internal quality control for quantitative measurements: principles and definitions, C24-A2. February 1999. National Committee for Clinical Laboratory Standards.
- 14 Krouwer, JS, Rabinowitz, R. How to improve estimates of imprecision. Clinical Chemistry, 1984, 30: 290-292.
- 15 Guideline-Measurement of Free Thyroid Hormones, C45A. 2004. National Committee for Clinical Laboratory Standards.

---

Beckman Coulter, Access, SYNCHRON LX, UniCel and DxI are trademarks of Beckman Coulter, Inc.; Beckman Coulter, Access, SYNCHRON LX, UniCel and DxI are registered in the USPTO and SIPO.

\*Lumi-Phos is a trademark of Lumigen, Inc, a subsidiary of Beckman Coulter, Inc.

\*\*Dynabeads is a registered trademark of Dynal A.S., Oslo, Norway.

\*\*\*ProClin is a trademark of Rohm and Haas Company or of its subsidiaries or affiliates.



Manufactured by:  
Beckman Coulter, Inc.  
250 S. Kraemer Blvd.  
Brea, CA 92821 U.S.A.



Beckman Coulter Ireland Inc.  
Mervue Business Park,  
Mervue, Galway,  
Ireland 353 91 774068

Printed in U.S.A.

Made in U.S.A. or Ireland and France

Revised July 2010

