Intended Use The Access Rubella IgG assay is a paramagnetic particle, chemiluminescent immunoassay for the qualitative and quantitative determination of IgG antibodies to the rubella virus in human serum using the Access Immunoassay Systems. The Access Rubella IgG assay aids in the diagnosis of rubella infection and the determination of immunity.

Summary and Explanation Rubella is a viral illness with worldwide distribution. The infection is usually benign or even unapparent in children or adults. Clinical manifestations include a generalized skin rash over the entire body, a low-grade fever, headache, and sometimes a sore throat. Infections in utero, particularly during the first four months of pregnancy, can lead to congenital defects such as deafness, cardiac problems, cataracts or glaucoma, and sometimes fetal death.^{1,2}

The detection of specific rubella antibodies is of great interest due to the teratogenic risk related to a rubella primary infection at the beginning of the pregnancy. Early methodologies used for antibody detection included serum neutralization, complement binding or immunofluorescence. These tests are difficult to perform and presented inherent problems of reproducibility. Subsequently, hemagglutination inhibition techniques allowed a more rapid diagnosis of both the acute infected state and patient immune status.^{3,4,5}

In 1971, Engvall and Perlmann described the first enzyme immunoassay procedures. The development of such methods have contributed to an improved specificity and sensitivity for research techniques with antigens and antibodies, including the rubella diagnosis field.⁵

Demonstration of rubella IgG antibody in a pregnant women prior to conception provides assurance of fetal protection from possible rubella viral infection during pregnancy. Vaccination efficiency is demonstrated by detection of rubella IgG antibody in serum following immunization. Appearance or significant rise of specific IgG concentration in two serum samples collected at least two weeks apart is indicative of rubella infection, even when typical symptoms may not be present.⁵

Principles of the Procedure The Access Rubella IgG assay is a enzyme immunoassay using an indirect technique. When a sample is added to a reaction vessel with paramagnetic particles coated with rubella viral membrane antigen,^{6,7,8,9} rubella antibodies bind to the antigen. After incubation in a reaction vessel, separation in a magnetic field and washing remove materials not bound to the solid phase. Alkaline phosphatase-conjugated monoclonal anti-human IgG antibody is then added and attaches to the IgG antibodies captured on the particles. A second separation and wash step removes unbound conjugate. A chemiluminescent substrate, Lumi-Phos* 530, is added to the reaction vessel and light generated by the reaction is measured with a luminometer. The light production is directly proportional to the concentration of IgG antibodies to the rubella virus in the sample. The amount of analyte in the sample is determined from a stored, multi-point calibration curve standardized against the WHO 2nd International Standard Preparation for Anti-Rubella Serum.

Product Information **Access Rubella IgG Reagent Pack**
Cat. No. 34430: 100 determinations, 2 packs, 50 tests/pack

- Provided ready to use.
- Store upright and refrigerate at 2 to 10°C.
- Refrigerate at 2 to 10°C for a minimum of two hours before use on the instrument.
- Stable until the expiration date stated on the label when stored at 2 to 10°C.

- Stable at 2 to 10°C for 28 days after initial use.
- Signs of possible deterioration are a broken elastomeric layer on the pack or control values out of range.
- If the reagent pack is damaged (i.e., broken elastomer), discard the pack.

R1a:	Paramagnetic particles coated with rubella (strain HPV 77) sucrose gradient purified antigen suspended in TRIS buffered saline, with bovine serum albumin (BSA), < 0.1% sodium azide, and 0.1% ProClin** 300.
R1b:	TRIS buffered saline with surfactant, BSA, < 0.1% sodium azide, and 0.1% ProClin 300.
R1c:	Mouse monoclonal anti-human IgG antibody (clone 125 A 15) - alkaline phosphatase (bovine) conjugate in TRIS buffered saline with surfactant, glycerol, BSA, mouse proteins, < 0.1% sodium azide.

Warnings and Precautions

- For *in vitro* diagnostic use.
- Patient samples and blood-derived products may be routinely processed with minimum risk using the procedure described. However, handle these products as potentially infectious according to universal precautions and good clinical laboratory practices, regardless of their origin, treatment, or prior certification. Use an appropriate disinfectant for decontamination. Store and dispose of these materials and their containers in accordance with local regulations and guidelines.
- Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal of liquids, flush with a large volume of water to prevent azide build-up.¹⁰
- Xi. Irritant: 0.1% ProClin 300.

R 43: May cause sensitization by skin contact.

S 28-37: After contact with skin, wash immediately with plenty of soap and water. Wear suitable gloves.

- The rubella antigen coated on the paramagnetic particles has been through a purification process. However, handle and dispose of reagent packs as potentially infectious material.
- The concentration of anti-rubella IgG in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity.
- The Material Safety Data Sheet (MSDS) is available upon request.

Specimen Collection and Preparation

1. Serum is the recommended sample. For the evaluation of immune status, collect a single serum sample. For the evaluation of seroconversion due to recent infection (the conversion of an individual patient's serum from negative to positive), collect two serum samples three weeks apart and test both sera within the same run.
2. Observe the following recommendations for handling, processing, and storing blood samples:¹¹
 - Collect all blood samples observing routine precautions for venipuncture.
 - Allow serum samples to clot completely before centrifugation.
 - Keep tubes stoppered at all times.
 - Within two hours after centrifugation, transfer at least 500 µL of cell-free sample to a tightly stoppered storage tube.
 - Store samples tightly stoppered at room temperature (15 to 30°C) for no longer than eight hours.
 - If the assay will not be completed within eight hours, refrigerate the samples at 2 to 8°C.
 - If the assay will not be completed within 48 hours, or for shipment of samples, freeze at -20°C or colder.
3. Use the following guidelines when preparing specimens:
 - Ensure residual fibrin and cellular matter have been removed prior to analysis.
 - Follow blood collection tube manufacturer's recommendations for centrifugation.

- Each laboratory should determine the acceptability of its own blood collection tubes and serum separation products. Variations in these products may exist between manufacturers and, at times, from lot-to-lot.
- Avoid repeated freezing and thawing of samples. A study of 10 serum samples ranging from 0–73 IU/mL (with six samples near the cutoff) indicates that samples may be frozen and thawed up to three times.

Materials Provided

R1 Access Rubella IgG Reagent Packs

Materials Required But Not Provided

- Access Rubella IgG Calibrators
Provided at zero, and approximately 10, 25, 50, 200 and 500 IU/mL.
Cat. No. 34435
- Access Rubella IgG Quality Control (QC) or other commercially available control material.
Provided as one negative and one low-level positive control (target mean of 22–43 IU/mL).
Cat. No. 34439
- Access Substrate
Cat. No. 81906
- Access, Access 2, SYNCHRON LXi:**
Access Wash Buffer II, Cat. No. A16792
UniCel DxI:
UniCel DxI Wash Buffer II, Cat. No. A16793

Procedural Comments

- Refer to the appropriate system manuals and/or Help system for a specific description of installation, start-up, principles of operation, system performance characteristics, operating instructions, calibration procedures, operational limitations and precautions, hazards, maintenance, and troubleshooting.
- Mix contents of new (unpunctured) reagent packs by gently inverting pack several times before loading on the instrument. Do not invert open (punctured) packs.
- Use twenty (20) µL of sample for each determination in addition to the sample container and system dead volumes. Refer to the appropriate system manuals and/or Help system for the minimum sample volume required.
- Report sample results in IU/mL or as reactive or non-reactive. See the “Results” section for further information on reporting results.

Procedure

Refer to the appropriate system manuals and/or Help system for information on managing samples, configuring tests, requesting tests, and reviewing test results.

Calibration Details

An active calibration curve is required for all tests. For the Access Rubella IgG assay, calibration is required every 28 days. Refer to the appropriate system manuals and/or Help system for information on calibration theory, configuring calibrators, calibrator test request entry, and reviewing calibration data.

Quality Control

Quality control materials simulate the characteristics of patient samples and are essential for monitoring the system performance of immunochemical assays. Because samples can be processed at any time in a “random access” format rather than a “batch” format, quality control materials should be included in each 24-hour time period.¹² More frequent use of controls or the use of additional controls is left to the discretion of the user based on good laboratory practices or laboratory accreditation requirements and applicable laws. Include Access Rubella IgG QC or other commercially available quality control materials that cover at least two levels of analyte. Follow manufacturer's instructions for reconstitution and storage. Include any additional controls if required for the laboratory's quality control program as guided by the appropriate regulatory agencies. Use controls in accordance with the appropriate accrediting organizations' requirements (for appropriate definitions of QC in the U.S.A., refer to CLSI

documents C24-A3, I/LA18-A2, and I/LA6-A). Users should refer to the appropriate system manuals and/or Help system for instructions on the use of the Quality Control functions and for selection of QC rules. In the U.S.A., it is suggested that a 1–3s QC rule be used for a low level reactive control. Quality control results that do not fall within acceptable ranges may indicate invalid test results. Examine all test results generated since obtaining the last acceptable quality control test point for this analyte.

Results Patient test results are determined automatically by the system software using a smoothing spline math model. The amount of analyte in the sample is determined from the measured light production by means of the stored calibration data. Patient test results can be reviewed using the appropriate screen. Refer to the appropriate system manuals and/or Help system for complete instructions on reviewing sample results.

The 15 IU/mL cutoff was based upon calibration to the WHO Second International Standard Preparation and a comparison study of 212 samples characterized by HAI. Because detection of nonimmune status is clinically more important, a 15 IU/mL cutoff was selected to assure the necessary specificity for the determination of immunity. Receiver-Operator Characteristic (ROC) curves supported the selection of the 15 IU/mL cutoff.

- All test samples < 10 IU/mL are considered non-reactive for the presence of rubella IgG antibodies. Patients with these results are considered to have an absence of immunity.
- Concentrations between ≥ 10 IU/mL and < 15 IU/mL are considered equivocal for determining immunity to rubella. Studies suggest that vaccinated individuals having these low levels of anti-rubella IgG do show a secondary immune response following re-vaccination but have not been challenged with wild rubella virus.⁵ A follow-up sample should be taken to further evaluate immune status. If the repeat sample is still equivocal, the sample may require testing by alternate methods.
- All test samples ≥ 15 IU/mL are considered to be reactive for the presence of rubella IgG antibodies and indicate acute or past infection, or vaccination.

The following is a recommended method for reporting the results obtained:

“The following results were obtained with the Access Rubella IgG EIA. Despite calibration by means of a reference preparation, values obtained with different manufacturer’s assay methods may not be used interchangeably. The magnitude of the reported IgG level cannot be correlated to an endpoint titer.”

Samples can be accurately measured within the linear range of the assay (approximately 10–500 IU/mL).

- If a sample contains < 10 IU/mL, report the results as non-reactive.
- If a sample contains ≥ 15 IU/mL, report as positive (or reactive) for anti-rubella IgG; presumed immune to rubella infection.
- If a sample contains more than the stated value of the highest Access Rubella IgG Calibrator (S5), report the result as greater than that value (i.e. approximately > 500 IU/mL). Alternatively, dilute one volume of sample with 9 volumes of Access Rubella IgG Calibrator S0 (zero). Refer to the appropriate system manuals and/or Help system for instructions on entering a sample dilution in a test request. The system reports the results adjusted for the dilution.

The meaning of IU/mL values above the cutoff and what constitutes a significant antibody level increase between acute and convalescent samples has not been determined.

For those patients suspected of early or recent seroconversion, a second serum sample should be collected three weeks later and evaluated along with the first serum sample for IgG antibody titer to rubella virus. For paired or serial serum samples, a conversion from a non-reactive to a reactive concentration of IgG antibody to rubella virus between the first and second serum samples should be considered as evidence of seroconversion due to recent infection. The laboratory should also evaluate the specimens for the presence of a significant level of IgM antibody to rubella to obtain additional serological data to aid in the diagnosis of a recent infection.

During the serological follow-up of rubella IgG non-reactive pregnant women, it is also recommended to measure anti-rubella IgM, as it is possible that anti-rubella IgG appearance may be slightly delayed with respect to that of anti-rubella IgM during or after recent infection.

Table 1: Examples of Single Serum Rubella IgG Test Results.

Serum Sample	Result (IU/mL)	Interpretation
Serum 1	0.4	Non-reactive; not Immune.
Serum 2	13.3	Equivocal: Rubella specific antibody present, evaluate further to determine immune status.
Serum 3	21.5	Reactive. Immune.
Serum 4	264.3	Reactive. Immune.
Serum 5	> 470.0	Reactive. Immune.

Table 2: Examples of Paired or Serial Serum Test Results.

Patient	Sample	Result (IU/mL)	Interpretation
1	1	0.1	Evidence of Seroconversion.
	2	310.1	
2	1	0.1	Evidence of Seroconversion.
	2	9.5	
	3	65.7	
	4	116.2	
3	1	0.2	No Evidence of Seroconversion.
	2	1.0	
4	1	55.6	No Evidence of Seroconversion.
	2	57.8	
5	1	0.4	Evidence of Seroconversion.
	2	58.9	
	3	118.9	

Limitations of the Procedure


1. The linear range of the assay is approximately 10–500 IU/mL. See the “Results” section for information on interpretation and reporting of results.
2. For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays. Additionally, other heterophile antibodies such as human anti-goat antibodies may be present in patient samples.^{13,14} Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.
3. The Access Rubella IgG results should be interpreted in light of the total clinical presentation of the patient, including: symptoms, clinical history, data from additional tests, and other appropriate information. The diagnosis of a recent infection by rubella virus can only be established on the basis of a combination of clinical and serological criteria. The result of a single serum sample does not constitute sufficient proof for the diagnosis of recent infection.
4. The presence of IgM antibody to rubella virus should also be evaluated as part of the serological surveillance of individuals suspected of rubella virus infection, as the appearance of rubella IgG antibody may occur slightly later than that of rubella IgM antibody.
5. Immunocompromised individuals and conditions, such as severe infection and immunosuppressive drug therapy, can result in the suppression of antibody levels below the detection threshold for the assay.
6. Assay performance characteristics have not been established for neonate or cord blood testing.
7. Assay performance characteristics have not been established using paired sera (acute and convalescent) to determine a significant change in IgG titer. For best results, pre-dilute paired sera with values greater than 200 IU/mL to values less than 200 IU/mL.
8. The meaning of IU/mL values above the cutoff and what constitutes a significant antibody increase between acute and convalescent samples has not been determined.
9. Assay performance characteristics have not been established for sera containing rheumatoid factor or ANA.

Expected Values

Rubella is worldwide in distribution and typically occurs more frequently in the spring and winter months. Incident rates vary with the epidemic cycle, the number of susceptibles among a population, and the interpersonal contact within the group. The disease is most prevalent among children of 5–9 years of age.² Serological studies in some areas show that up to 95% of pregnant women are seropositive. Several countries have implemented a routine immunization program for infants at the age of approximately 15 months.

The prevalence of reactive results for anti-rubella IgG antibodies was 85% as determined from 414 serum samples obtained from a normal Memphis, TN (U.S.A.) blood donor population using the Access Rubella IgG assay. The distribution of results is provided in the following histogram.

Access Rubella IgG: Screening of U.S. blood donor population.

Specific Performance Characteristics

Methods Comparison

Two independent sites in France evaluated the clinical performance of the Access Rubella IgG assay by comparison to HAI. A total of 784 patients were tested. For site 1, 399 of the samples were from pregnant women followed within the scope of the compulsory serological survey of pregnant women in France which includes monthly monitoring of those testing as seronegative. Eleven of the site 1 samples (all pregnant) were equivocal and removed from the results calculations. For site 2, 290 samples were from pregnant women. Ten of the site 2 samples (all female, 8 pregnant) were equivocal and removed from the results calculations. The prospective samples were from a normal population consisting of 28 pregnant women, 21 nonpregnant women and 1 man. Retrospective (frozen) samples were utilized in these studies to expedite obtaining the desired number of reactive and non-reactive samples to sufficiently evaluate the assay. These samples may have undergone multiple freeze/thaw cycles. In an internal study, the performance of the assay was found to be unaffected by fresh vs. frozen/thawed samples (up to 3 cycles). The relative sensitivities, specificities, and the corresponding 95% confidence intervals are provided in Table 3 for fresh and frozen samples. The stated estimates of sensitivity are the best estimates of the “relative” sensitivities of this product given the data collected. A 95% confidence interval constructed from such data can be expected to contain the “relative” sensitivity with a probability of 95%. The widths of the stated 95% confidence intervals of sensitivity almost exclusively depend on the number of positives sampled.

Table 3: Access versus HAI.


		HAI							Relative Sensitivity (%)	95% Conf. Int.	Relative Specificity (%)	95% Conf. Int.
n	Access	+	+	+	–	–	–					
Site 1	385	Frozen	280	0	1	10	1	93	100	98.7–100	99	94.2–100
Site 1	25	Fresh	22	0	0	0	0	3	100	84.6–100	100	29.2–100
Site 2	349	Frozen	264	10	10	0	0	65	96	93.4–98.2	100	94.5–100
Site 2	25	Fresh	20	0	0	0	0	5	100	83.2–100	100	47.8–100
TOTAL	784		586	10	11	10	1	166	98	96.9–99.2	99	96.7–100

Table 4 shows the results of the pregnant and non-pregnant populations (male and female) from both sites. At site 1, there was one newborn patient and one patient of unknown sex that were not included. These two samples were not discordant and were reactive by both Access and HAI.

Table 4: Pregnant/Non-Pregnant Population.

		HAI	+	+	+	-	-	-	Relative Sensitivity (%)	95% Conf. Int.	Relative Specificity (%)	95% Conf. Int.
	n	Access	+	-	EQ	EQ	+	-				
Site 1	399	Pregnant	291	0	1	10	1	96	100	98.7–100	99	94.4–100
Site 2	290	Pregnant	216	10	8	0	0	56	96	92–97.9	100	93.6–100
Site 1	5	Female	5	0	0	0	0	0	100	47.8–100	NA	NA
Site 1	4	Male	4	0	0	0	0	0	100	39.8–100	NA	NA
Site 2	75	Female	64	0	2	0	0	9	100	94.4–100	100	66.4–100
Site 2	9	Male	4	0	0	0	0	5	100	39.8–100	100	47.8–100
TOTAL	782		584	10	11	10	1	166	98	96.9–99.2	99	96.7–100

The distribution of all 784 clinical samples from the comparison to HAI is provided in the histogram below. There were 177 HAI negative and 607 HAI positive samples.

The Access Rubella IgG assay was compared to another commercially available EIA using 306 retrospective samples at a western U.S. site. There were 181 samples identified as reactive by both tests and 79 samples identified as non-reactive by both tests. The 46 remaining samples were EIA reactive /Access non-reactive. Of the 44 discordant samples further analyzed by a latex agglutination assay, 25 agreed with Access.

CDC Panel Results

The following information is from a serum panel obtained from the CDC and tested by an external site using the Access Rubella IgG assay. The results are presented as a means to convey further information on the performance of this assay with a masked, characterized serum panel. This does not imply an endorsement of the assay by the CDC.

The panel consists of 82 positive sera and 18 negative sera. The Access Rubella IgG assay demonstrated 100% total agreement with the CDC results, there was 100% agreement with the positive specimens and 100% agreement with the negative specimens.

CDC Biological Standard Results

The low titer (21.0 IU/mL) anti-rubella human reference serum CDC biological standard was tested neat and diluted 1:2 as described in the CLSI document I/LA6-A at three external sites. The average of the results from the three sites using the Access Rubella IgG assay were 25.1 IU/mL for the neat sample and 11.8 IU/mL for the 1:2 dilution.

Congenital Rubella

Serum samples related to clinically documented cases of congenital infection by the rubella virus are extremely rare, however, the following data obtained during the external evaluations of the assay demonstrate good agreement between clinical information and assay results obtained with the Access Rubella IgG assay.

Table 5: Congenital Rubella.

Subject	Draw Date	Clinical Information	Access Rubella IgG (IU/mL)	HAI Titer	Rubella IgM
Case 1					
Mother	8/13/93	10 days after rash	275	≥ 512	positive
Newborn	9/30/93	congenital rubella	> 470	≥ 512	positive
Case 2					
Mother	7/6/93	70 days after rash	375	≥ 512	equivocal
Newborn	9/15/93	congenital rubella	> 470	≥ 512	positive

Immune Response Studies

Multiple samples drawn from patients with a recent rubella infection or vaccination were tested to evaluate the immune response. Evidence of acute rubella infection was supported in 38 out of 44 of these individuals by a reactive Rubella IgM antibody test result. Site 1 tested 23 samples from 15 patients with recent rubella infection and 22 samples from 9 patients recently vaccinated. Site 2 tested 21 samples from 10 patients with recent rubella infection and 15 samples from 7 patients recently vaccinated. Site 3 tested 21 seroconversion samples from 3 patient panels purchased from serum suppliers. The increase in rubella IgG as measured by the Access Rubella IgG assay correlated to increasing titers of rubella IgG as detected by the comparative rubella IgG assay. Of the 74 reactive samples as determined by the comparison method, the Access Rubella IgG identified 67 as reactive, 1 as equivocal and 6 as non-reactive. Of the 28 non-reactive samples as determined by the comparison method, the Access Rubella IgG identified 2 as equivocal and 26 as non-reactive. Clinical information and further testing results available on the 9 discordant samples are shown in Table 6.

Table 6: Discordant Sample Summary.

Patient	Clinical Information	HAI Titer [†]	Access Rubella IgG EIA Result	Second EIA Rubella IgG Result	Rubella IgM Result	Latex Agglutination Result [†]
455	3 days after rash	1:32 positive	12.1 IU/mL equivocal	35 IU/mL positive	positive	> 10 positive
BOU	had rash	1:32 positive	non-reactive	not done	positive	> 10 positive
456	before rash	1:8 positive	non-reactive	not done	negative	< 10 negative
461	before rash	1:8 positive	non-reactive	4.5 IU/mL negative	negative	< 10 negative
450	not available	1:8 positive	non-reactive	not done	negative	> 10 positive
458	before rash	1:32 positive	non-reactive	0.8 IU/mL negative	equivocal	< 10 negative
447	before vaccination	1:128 positive	non-reactive	0.0 IU/mL negative	negative	< 10 negative
T37581	not available	not done	13.8 IU/mL equivocal	negative	positive	not done
G40955	not available	not done	13.6 IU/mL equivocal	negative	positive	not done

[†] HAI and Latex Agglutination detect both IgM and IgG

Samples BOU and 455 are from patients with rash or just three days after rash. They are IgM positive and too early in seroconversion to be reactive for IgG. As the Access result was non-reactive on the before rash sample for patient 455, the Access equivocal result on the three days after rash sample indicates the IgG titer is rising.

Samples 456, 461, 458, and 447 are from patients before rash or vaccination. No other test results on these samples confirm the HAI positive result.

Sample 450 is from very early seroconversion as indicated by the positive (at the cutoff of 1:8) HAI and latex agglutination results. The IgG level is rising (9.1 IU/mL) but not yet positive. The next bleed from

this individual two months later, was definitely positive having a 1:256 HAI titer and 210 IU/mL by Access Rubella IgG.


For samples T37581 and G40955, this bleed is the one just before the bleed in which the second EIA converts from negative to positive. The IgM results are positive indicating seroconversion is in progress.

Considering the very early seroconversions in this sample population and the fact that IgG levels appear after IgM levels which are detected by HAI, the performance of the Access Rubella IgG assay on immune response samples was very good.

Dilution Recovery (Linearity)

Linearity with WHO Standard


Five dilutions of the WHO second International Standard for Anti-Rubella Serum were tested in triplicate. One of the three replicates of the 400 IU/mL dilution was over range. In order to avoid biasing the linear regression, none of the replicates from this dilution are included in the generation of the line graph below.

Linearity with “Over Range” Patient Samples

Five “over range” patient samples were diluted 1:2, 1:4, and 1:8 and tested in triplicate in three different assays. Expected values for each patient sample were the means of the 1:2 dilution. All observed values for the five samples are presented below. The regression line is calculated using data from only those dilutions which had no “over range” sample results.

Access Rubella IgG: Five Samples and Dilutions

Imprecision

Three sites (one internal and two external) assayed two controls and four samples in triplicate, twice a day for three days. One calibration was performed at the beginning of the study. The imprecision of the non-reactive sample was calculated from the relative light units (RLU), all other calculations used the IU/mL values. Results are summarized in Table 7.

Table 7: Rubella IgG Three Site Imprecision.

Level	n	Mean	Within Run		Within Site		Total Imprecision	
			SD	% CV	SD	% CV	SD	% CV
QC1 (RLU)	54	44,139	6514	14.8	6426	14.6	7178	16.3
QC2	54	43.9	2.5	5.7	2.6	5.9	2.6	5.9
P2	53	286.5	26.0	9.1	30.7	10.7	33.3	11.6
P3	54	416.4	27.7	6.7	45.6	11.0	48.4	11.6
P4	54	16.8	0.6	3.8	1.0	6.0	1.4	8.2
P6	54	75.6	2.5	3.2	4.1	5.4	4.4	5.8

Lot to Lot Reproducibility

For each product lot below, one replicate of 20 non-reactive patient samples (selected from a bank of 60 non-reactive patient samples) and duplicates of five reactive patient samples were tested. All 60 non-reactive patient samples tested were found to be non-reactive in the Access Rubella IgG assay. The results for the reactive patient samples (in IU/mL) are presented below.

Table 8: Lot to Lot Reproducibility.

Lot	Patient Sample				
	P16	P2	P3	P9	P13
490211	14	52	119	280	> 459
	16	52	126	245	> 459
490213	13	47	123	225	> 458
	12	49	111	251	> 458
490214	13	45	106	195	450
	12	46	102	202	> 458
490215	15	60	129	283	>458
	15	59	138	306	> 458
590216	17	65	199	373	> 500
	14	56	156	364	> 500

Analytical Specificity/Interferences

Samples artificially spiked to contain up to 9 g/dL albumin, up to 20 mg/dL bilirubin, up to 3600 mg/dL triolein, and up to 2000 mg/dL hemoglobin do not affect the measured concentration of the detected rubella IgG antibodies.

The following samples were evaluated and found non-reactive in the Access Rubella IgG assay indicating a lack of cross-reactivity and/or non-specific reactivity with these specimens. The non-reactive Access Rubella IgG results were confirmed by testing with another rubella IgG EIA method.

Table 9: Cross Reactivity and Interference

Number	Sample type
16	CMV IgG samples
17	EBV-VCA IgG samples
2	EBV-VCA IgM samples
16	HSV 1 IgG samples
12	HSV 2 IgG samples
5	Measles IgG samples
7	Measles IgM samples
7	Toxo IgG samples

RUBELLA IgG CALIBRATORS

REF 34435

Intended Use The Access Rubella IgG Calibrators are intended to calibrate the Access Rubella IgG assay for the qualitative and quantitative determination of IgG antibodies to the rubella virus in human serum using the Access Immunoassay Systems.

Summary and Explanation Quantitative assay calibration is the process by which samples with known analyte concentrations (i.e. assay calibrators) are tested like patient samples to measure the response. The mathematical relationship between the measured responses and the known analyte concentrations establishes the calibration curve. This mathematical relationship, or calibration curve, is used to convert RLU (Relative Light Unit) measurements of patient samples to specific quantitative analyte concentrations.

Traceability The measurand (analyte) in the Access Rubella IgG Calibrators is traceable to the WHO Second International Standard Preparation for Anti-Rubella Serum (2nd ISP). Traceability process is based on EN ISO 17511.

The assigned values were established using representative samples from this lot of calibrator and are specific to the assay methodologies of the Access reagents. Values assigned by other methodologies may be different. Such differences, if present, may be caused by inter-method bias.

Product Information **Access Rubella IgG Calibrators**
Cat. No. 34435: S0–S5, 1.0 mL/vial

- Provided ready to use.
- Store upright and refrigerate at 2 to 10°C.
- Mix contents by gently inverting before use. Avoid bubble formation.
- Stable until the expiration date stated on the label when stored at 2 to 10°C.
- Open vial stability is typically until the expiration date stated on the vial label when properly stored and handled.
- Signs of possible deterioration are control values out of range.
- Refer to calibration card for exact concentrations.

S0:	Equine serum with 0 IU/mL anti-rubella IgG, and < 0.1% sodium azide.
S1, S2, S3, S4, S5:	Equine serum and human defibrinated plasma containing approximately 10, 25, 50, 200 and 500 IU/mL of human anti-rubella IgG and < 0.1% sodium azide.
Calibration Card:	1

Warnings and Precautions

- For *in vitro* diagnostic use.
- Patient samples and blood-derived products may be routinely processed with minimum risk using the procedure described. However, handle these products as potentially infectious according to universal precautions and good clinical laboratory practices, regardless of their origin, treatment, or prior certification. Use an appropriate disinfectant for decontamination. Store and dispose of these materials and their containers in accordance with local regulations and guidelines.

- Human source material used in the preparation of the reagent has been tested and found negative or non-reactive for Hepatitis B, Hepatitis C (HCV), and Human Immunodeficiency Virus (HIV-1 and HIV-2). Because no known test method can offer complete assurance that infectious agents are absent, handle reagents and patient samples as if capable of transmitting infectious disease.¹⁵
 - Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal of liquids, flush with a large volume of water to prevent azide build-up.¹⁰
 - The Material Safety Data Sheet (MSDS) is available upon request.
-

Procedure Refer to the appropriate system manuals and/or Help system for information on calibration theory, configuring calibrators, calibrator test request entry, and reviewing calibration data.

Calibration Details The Access Rubella IgG Calibrators are provided at six levels - zero and approximately 10, 25, 50, 200 and 500 IU/mL - prepared from equine serum and human defibrinated plasma positive for anti-rubella IgG. Assay calibration data are valid up to 28 days.
Calibrators run in duplicate.

Limitations of the Procedure If there is evidence of microbial contamination or excessive turbidity in a reagent, discard the vial.

RUBELLA IgG QC

REF 34439

Intended Use The Access Rubella IgG QC is intended for monitoring system performance of the Access Rubella IgG assay.

Summary and Explanation Quality control materials simulate the characteristics of patient samples and are essential for monitoring the system performance of the Access Rubella IgG immunoassay. In addition, they are an integral part of good laboratory practices.^{12,16,17,18,19,20,21} When performing assays with Access reagents for IgG antibodies to the rubella virus, include quality control materials to validate the integrity of the assays. The assayed values should fall within the acceptable range if the test system is working properly.

Traceability The measurand (analyte) in the Access Rubella IgG QC is traceable to the WHO Second International Standard Preparation for Anti-Rubella Serum (2nd ISP). Traceability process is based on EN ISO 17511. The assigned values were established using representative samples from this lot of QC and are specific to the assay methodologies of the Access reagents. Values assigned by other methodologies may be different. Such differences, if present, may be caused by inter-method bias.

Product Information **Access Rubella IgG QC**
Cat. No. 34439: 2.5 mL/vial, 3 vials each level

- Provided ready to use.
- Store upright and refrigerate at 2 to 10°C.
- Mix contents by gently inverting before use. Avoid bubble formation.
- Stable until the expiration date stated on the label when stored at 2 to 10°C.
- After initial use, vials are stable for 30 days when correctly handled and stored.
- Signs of possible deterioration are control values out of range.
- Refer to the QC value card for mean values and standard deviations (SD).

QC 1:	Human defibrinated plasma with < 0.1% sodium azide; contains no detectable level of anti-rubella IgG as assayed using the Access Rubella IgG assay.
QC 2:	Human defibrinated plasma with < 0.1% sodium azide; contains a low level of anti-rubella IgG (target mean of 22–43 IU/mL) as assayed using the Access Rubella IgG assay.
QC Value Card:	1

Warnings and Precautions

- For *in vitro* diagnostic use.
- Patient samples and blood-derived products may be routinely processed with minimum risk using the procedure described. However, handle these products as potentially infectious according to universal precautions and good clinical laboratory practices, regardless of their origin, treatment, or prior certification. Use an appropriate disinfectant for decontamination. Store and dispose of these materials and their containers in accordance with local regulations and guidelines.
- Human source material used in the preparation of the reagent has been tested and found negative or non-reactive for Hepatitis B, Hepatitis C (HCV), and Human Immunodeficiency Virus (HIV-1 and

HIV-2). Because no known test method can offer complete assurance that infectious agents are absent, handle reagents and patient samples as if capable of transmitting infectious disease.¹⁵

- Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal of liquids, flush with a large volume of water to prevent azide build-up.¹⁰
 - The Material Safety Data Sheet (MSDS) is available upon request.
-

Procedure Determine the concentration of IgG antibodies to the rubella virus in the Access Rubella IgG QC materials using the Access Immunoassay System in the same manner as a patient sample. Because samples can be processed at any time in a “random access” format rather than a “batch” format, quality control materials should be included in each 24-hour time period.¹² More frequent use of controls or the use of additional controls is left to the discretion of the user based on good laboratory practices or laboratory accreditation requirements and applicable laws. Refer to the appropriate system manuals and/or Help system for information on quality control theory, configuring controls, quality control sample test request entry, and reviewing quality control data.

Limitations of the Procedure

1. Quality controls must be run within the 24 hour period prior to running patient samples. Include any additional controls if required for the laboratory's quality control program as guided by the appropriate regulatory agencies.
2. Use quality controls in accordance with the appropriate accrediting organizations' requirements (for appropriate definitions of QC in the U.S.A., refer to CLSI documents C24-A3, I/LA18-A2, and I/LA6-A). Users should refer to the appropriate system manuals and/or Help system for instructions on the use of the Quality Control functions and for selection of QC rules.
3. In the U.S.A., it is suggested that a 1–3s QC rule be used for a low level reactive control.
4. Quality control results that do not fall within acceptable ranges may indicate invalid test results. Examine all test results generated since obtaining the last acceptable quality control test point for this analyte.
5. If there is evidence of microbial contamination or excessive turbidity in a reagent, discard the vial.

Expected Values The expected means (**0**) and SD's (σ) for the Access Rubella IgG QC controls (QC1 and QC2) are provided on the QC value card contained in the kit for initial Access Rubella IgG assay quality control system configuration. Each laboratory should establish its own acceptability criteria by selecting the QC rules to be applied to the control results. Individual control results should fall within the initial acceptable range, however, each laboratory should update the mean and SD after sufficient data has been collected.^{20,21}

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Code 34430: Coffret de 100 dosages, 2 packs, 50 tests/pack

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- 1 Wolinsky JS. Rubella. In Virology, 2nd Ed. 815-38. Edited by Fields BN, et al. New York. Raven Press, Ltd., 1990.
- 2 Assaad F, Ljungars-Esteves K. Rubella – world impact. Rev Infect Dis 1985; Vol. 7 (Sup 1): S29-S36.
- 3 Boué A, Nicolas A, Lang R. La réaction d'inhibition de l'hémagglutination pour la sérologie de la rubéole. Ann Inst Pasteur 1968; 114: 317-30.
- 4 Forsgren, M. Standardization of techniques and reagents for the study of rubella antibody. Rev Infect Dis 1985; 7 (Suppl. 1): 129-32.
- 5 Approved Guideline – Detection and Quantitation of Rubella IgG Antibody: Evaluation and Performance Criteria for Multiple Component Test Products, Specimen Handling, and Use of Test Products in the Clinical Laboratory, I/LA6-A. 1997. Clinical and Laboratory Standards Institute.
- 6 Dorsett PH, Miller DC, Green K, Byrd RF. Structure and function of the rubella virus proteins. Rev Infect Dis 1985; 7 (Suppl. 1): 150-56.
- 7 Kalkkinen N, Oker-Blom C, Pettersson RF. Three genes code for rubella virus structural proteins E1, E2a, E2b. J Gen Virol 1984; 65: 1549-57.
- 8 Pettersson R, et al. Molecular and antigenic characteristics and synthesis of rubella virus structural proteins. Rev Infect Dis 1985; 7 (Suppl. 1): 140-49.
- 9 Lucas G, et al. Serological diagnosis of anti-rubella immunoglobulines (IgG) by an immunoenzymatic assay, commercially available as a kit. Nice: 4th European Congress of Clinical Microbiology 1989; 308: 20.
- 10 DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard. Available <http://www.cdc.gov/niosh>.
- 11 Approved Guideline – Procedures for the Handling and Processing of Blood Specimens, H18-A3. 2004. Clinical and Laboratory Standards Institute.
- 12 Cembrowski GS, Carey RN. Laboratory quality management: QC=QA. ASCP Press, Chicago, IL, 1989.
- 13 Kricka L. Interferences in immunoassays – still a threat. Clin Chem 2000; 46: 1037–1038.
- 14 Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613–621.
- 15 HHS Publication, 4th ed., April 1999. Biosafety in Microbiological and Biomedical Laboratories. Available <http://www.cdc.gov/od/ohs/biosfty/bmb14/bmb14toc.htm>
- 16 Broome HE, Cembrowski GS, Kahn SN, Martin PL, Patrick CA. Implementation and use of a manual multi-rule quality control procedure. Lab Med 1985; 16: 533-537.
- 17 Westgard JO, Barry PL, Hunt MR, Groth T. A multi-rule Shewhart chart for quality control in clinical chemistry. Clin Chem 1981; 27: 493-501.
- 18 Koch DD, Oryall JJ, Quam EF, Feldbruegger DH, et al. Selection of medically useful QC procedures for individual tests done in a multitest analytical system. Clin Chem 1990; 36: 230-233.
- 19 Mugan K, Carlson IH, Westgard JO. Planning QC procedures for immunoassays. J Clin Immunoassay 1994; 17: 216-222.
- 20 Approved Guideline – Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions, C24-A3. June 2006. Clinical and Laboratory Standards Institute.
- 21 Garrett PE. Quality is quantitative: so how do we QC qualitative tests? J Clin Immunoassay 1994; Vol 17, No. 4: 231.

Beckman Coulter, the stylized logo, Access, SYNCHRON LX, UniCel and DxI are trademarks of Beckman Coulter, Inc. and are registered in the USPTO.

*Lumi-Phos is a trademark of Lumigen, Inc., a subsidiary of Beckman Coulter, Inc./ *Lumi-Phos est une marque déposée de Lumigen, Inc., une filiale de Beckman Coulter, Inc./ *Lumi-Phos es una marca registrada de Lumigen, Inc., una subsidiaria de Beckman Coulter, Inc./ *Lumi-Phos è un marchio di Lumigen, Inc., un'affiliata di Beckman Coulter, Inc./ *Lumi-Phos ist eine Marke von Lumigen, Inc., einer Tochtergesellschaft von Beckman Coulter, Inc./ *Lumi-Phos é uma marca registrada da Lumigen, Inc., uma companhia subsidiária da Beckman Coulter, Inc./ *Lumi-Phos 는 Beckman Coulter, Inc. 의 자회사인 Lumigen, Inc. 의 상표입니다 .

**ProClin is a trademark of Rohm and Haas Company or of its subsidiaries or affiliates./ **ProClin est une marque déposée de Rohm and Haas Company ou de ses filiales./ **ProClin es una marca registrada de Rohm and Haas Company o sus filiales o asociadas./ **ProClin è un marchio di Rohm and Haas Company, o delle sue filiali o controllate./ **ProClin ist eine Marke der Rohm and Haas Company bzw. deren Tochtergesellschaften oder Partnerfirmen./ **ProClin é uma marca registrada da Rohm and Haas Company ou de suas subsidiárias e filiais./ **ProClin 은 Rohm and Haas Company 또는 자회사나 계열사의 상표입니다 .

Manufactured for:/Fabriqué pour:/Fabricado para:/ Prodotto
per:/Hergestellt für:/Fabricado para:/ 제조업체 :
Beckman Coulter, Inc.
250 S. Kraemer Blvd.
Brea, CA 92821 U.S.A.

Beckman Coulter do Brasil Com e Imp de Prod Lab Ltda
Estr dos Romeiros, 220 - Galpao G3 - KM 38.5
06501-001 - Sao Paulo - SP - Brasil
CNPJ: 42.160.812/0001-44

EC REP

Beckman Coulter Ireland Inc.
Mervue Business Park,
Mervue, Galway,
Ireland 353 91 774068

Printed in U.S.A./Imprimé aux Etats-Unis/Impreso en U.S.A./ Stampato negli Stati Uniti d'America/Gedruckt in den USA./Impresso nos Estados
Unidos da América./ 미국에서 인쇄

Made in France/Fabriqué en France/Fabricado en Francia/Fabbricato in Francia/Hergestellt in Frankreich/Fabricado na França./
프랑스에서 제작

Revised November 2013/Révision en novembre 2013/Revisado en noviembre de 2013/Revisione novembre 2013/Überarbeitet November 2013/
Revisado em Novembro 2013/2013 년 11 월 개정