

## SUBSTRATE

**REF 81906**

**Intended Use** Access Substrate is intended for use with the Access Immunoassay Systems and specific Access immunoassay reagents.

**Product Information**

- Provided ready to use.
- Refer to the following chart for storage conditions and stability. An increase in substrate background measurements may indicate instability.

Condition	Storage	Stability
Unopened	2–8°C	Until expiration date stated on the label
Equilibration prior to use (unopened)	15–30°C (room temperature)	Minimum 18 hours Maximum 14 days
In use (opened)	External fluids tray substrate position	Maximum 14 days

**R2 Access Substrate: 4 x 130 mL.**

Lumi-Phos\* 530 (buffered solution containing dioxetane Lumigen\* PPD, fluorescer, and surfactant).

**Warnings and Precautions**

- For *in vitro* diagnostic use.
- Substrate is sensitive to air exposure. Keep tightly closed at all times. Do not pool bottles of substrate.

**Directions for Use** Refer to the appropriate system Manuals and/or Help system for detailed instructions.

Beckman Coulter and Access are trademarks of Beckman Coulter, Inc.; Beckman Coulter and Access are registered in the USPTO and SIPO.

\*Lumi-Phos and Lumigen are trademarks of Lumigen, Inc., a subsidiary of Beckman Coulter, Inc.



Manufactured by:  
Beckman Coulter, Inc.  
250 S. Kraemer Blvd.  
Brea, CA 92821 U.S.A.

Printed in U.S.A.  
Made in U.S.A.  
Revised June 2010

**EC REP**

Beckman Coulter Ireland Inc.  
Mervue Business Park,  
Mervue, Galway,  
Ireland 353 91 774068

