

SYSTEM CHECK SOLUTION

REF 81910

Intended Use For use with the Access Immunoassay Systems in the weekly maintenance System Check Procedure.

Product Information

- Provided ready to use.
- Stable until the expiration date stated on the label when stored at 2 to 8°C.

R4 Access System Check Solution: 6 x 4.0 mL.
Alkaline phosphatase, 1% bovine serum albumin (BSA), 0.25% ProClin* 300, < 0.1% sodium azide.

Warnings and Precautions

- For *in vitro* diagnostic use.
- Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal of liquids, flush with a large volume of water to prevent azide build-up.¹
- Xi. Irritant: 0.25% ProClin 300.



R 43: May cause sensitization by skin contact.

S 28-37: After contact with skin, wash immediately with plenty of soap and water. Wear suitable gloves.

Directions for Use Refer to the appropriate system Manuals and/or Help system for detailed instructions.

References

- 1 DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard. Available <http://www.cdc.gov/niosh>.

Beckman Coulter and Access are trademarks of Beckman Coulter, Inc.; Beckman Coulter and Access are registered in the USPTO and SIPO.

*ProClin is a trademark of Rohm and Haas Company or of its subsidiaries or affiliates.



Manufactured by:
Beckman Coulter, Inc.
250 S. Kraemer Blvd.
Brea, CA 92821 U.S.A.

Printed in U.S.A.
Made in U.S.A.
Revised July 2010



Beckman Coulter Ireland Inc.
Mervue Business Park,
Mervue, Galway,
Ireland 353 91 774068

